

# TUBE STRIPPER

## Product Features

High quality aluminum alloy  
Anti-rust, light weight and sturdy  
Ease of cleaning by wiping or soaking in IPA



Ergonomic design for good grip and stripping  
Ease fatigue over prolong use




Self-centering roller system keeps tubing in position  
Smooth traction for effective stripping of blood components



Internal spring system keeps handles in open position  
Minimize entanglement with blood tubing



<b>Model</b>	<b>HS - 002</b>
Base Material	Aluminum Alloy & SUS
Rollers	Special Composite Steel, Non-Rust
Cutters	Not Available
Features	Self-Centering Rollers
Cutter	None
Net Weight (g) +/- 5%	68



**HS-002**