

TERUFLEX® and IMUFLEX® Blood Bags

All blood bag systems feature an XT150 primary container. Double systems have one 400mL satellite bag. TERUFLEX triple systems have one 400mL satellite bag and one XT-612® platelet storage bag. Quadruple systems have two 400mL satellite bags and one XT-612 platelet storage bag. Product features include:

- 16 gauge ultra-thinwall needle for donor comfort
- Kink resistant tubing
- CLIKTIP® effortless inline closure device
- Easy to open outlet port peel tabs
- Latex free
- Color coded blister packaging

Legend for Blood Bag Sampling Arms:

A6	IMUFLEX Blood Bags with Diversion Blood Sampling Arm®, DonorCare ® Needle Guard and pre-attached luer adapter and tube holder. Features an integral sampling pouch with pre-attached luer adapter and tube holder to collect blood for test samples from the donor prior to collecting the unit. These bags feature a DonorCare Needle Guard pre-attached.
A2	TERUFLEX and IMUFLEX Blood Bags with Diversion Blood Sampling Arm and DonorCare Needle Guard. Features an integral sampling pouch to collect blood for test samples from the donor prior to collecting the unit. These bags feature a DonorCare Needle Guard pre-attached.
A3	TERUFLEX Blood Bags with Blood Sampling Arm and DonorCare Needle Guard. The integral sampling arm system is easy to use for obtaining donor samples after the unit is collected. Also features the DonorCare Needle Guard pre-attached.
A4	TERUFLEX Blood Bags with Blood Sampling Arm. The integral sampling arm system is easy to use for obtaining donor samples after the unit is collected.
A	TERUFLEX Blood Bags with Standard Donor Tubing. Straight line tubing, without a sampling system.

IMUFLEX WB-SP Blood Bag System with in-line Leukocyte Reduction Filter with Diversion Blood Sampling Arm and pre-attached luer adapter and tube holder

For the preparation of leukocytes reduced Red Blood Cells, Platelets, and Plasma from Whole Blood. Contains a single in-line **IMUFLEX** leukocyte reduction filter connected to a 500mL XT150 primary container with one 400mL satellite bag and one XT-612 Platelet bag for component storage. Donor tubing has a pre-attached **DonorCare** Needle Guard.

CPD/OPTISOL® Anticoagulant and Additive Solution

500 mL Collection Bags **Catalog No.**
 Triple Blood Bag **1BB*LGQ506A6** 15/case



IMUFLEX WB-RP Blood Bag System with in-line Leukocyte Reduction Filter and Diversion Blood Sampling Arm

For the preparation of leukocyte reduced whole blood. Contains an in-line **IMUFLEX** leukocyte reduction filter connected to a 500mL XT150 primary container with two 400mL satellite bags for component storage. Not for the preparation of platelets. Donor tubing has a pre-attached **DonorCare** Needle Guard.

CPD/OPTISOL Anticoagulant and Additive Solution

500 mL Collection Bags **Catalog No.**
 Triple Blood Bag **1BB*WGQ506A2** 18/case



- A6-** Diversion Blood Sampling Arm, DonorCare Needle Guard and pre-attached tube holder and luer adapter
- A2-** Diversion Blood Sampling Arm and DonorCare Needle Guard

TERUFLEX® Blood Bag Systems

CPD/OPTISOL® Anticoagulant and Additive Solution

450 mL Collection Bags	Catalog No.	
Double Blood Bag	1BB*AGD456A	24/case
	1BB*AGD456A2	24/case
	1BB*AGD456A3	24/case
Triple Blood Bag	1BB*AGT456A	24/case
	1BB*AGT456A2	24/case
	1BB*AGT456A3	24/case
Quadruple Blood Bag	1BB*AGQ456A	24/case
	1BB*AGQ456A2	24/case
500 mL Collection Bags	Catalog No.	
Double Blood Bag	1BB*AGD506A2	24/case
	1BB*AGD506A3	24/case
Triple Blood Bag	1BB*AGT506A2	24/case
	1BB*AGT506A3	24/case
Quadruple Blood Bag	1BB*AGQ506A2	18/case



CPDA-1 Anticoagulant

450 mL Collection Bags	Catalog No.	
Single Blood Bag	1BB*SCD456A	60/case
	1BB*SCD456A3	36/case
	1BB*SCD456A4	48/case
Double Blood Bag	1BB*DCD456A	30/case
	1BB*DCD456A3	30/case
	1BB*DCD456A4	30/case
Triple Blood Bag	1BB*TCD456A4	30/case
Quadruple Blood Bag	1BB*QCD456A3	24/case
500 mL Collection Bags	Catalog No.	
Single Blood Bag	1BB*SCD506A3	36/case
Double Blood Bag	1BB*DCD506A3	30/case
Quadruple Blood Bag	1BB*QCD506A3	24/case



- A2- Diversion Blood Sampling Arm® and DonorCare® Needle Guard
- A3- Blood Sampling Arm and DonorCare Needle Guard
- A4- Blood Sampling Arm
- A- No Sampling System


TERUFLEX® Dry Blood Bag

Product Description	Catalog No.	
Empty blood bag with 16 gauge ultra-thinwall donor needle, 600mL XT150 container.	1BB*D606A	100/case



TERUFLEX Transfer Bags

Product Description	Catalog No.	
150mL, One Coupler	1BB*T015CB70	100/case
300mL, One Coupler	1BB*T030CB71	100/case
600mL, One Coupler	1BB*T060CB71	100/case
800mL, Two Couplers	1BB*T080BB71	100/case
1000mL, Two Couplers	1BB*T100BB71	100/case
2000mL, Two Couplers	1BB*T200BB71	100/case



TERUFUSION® Blood Transfusion Sets/Miscellaneous (Latex Free)

Product Description	Catalog No.	
Y-type blood administration set with a 170 micron nylon mesh filter, and needle/catheter adapter, approximate length 60"	1TY*A2005020	50/case
Straight-type blood administration set with a 170 micron nylon mesh filter, and needle/catheter adapter, approximate length 58". Includes one 18 ga. x 1.5" needle	1TB*A200B	50/case
Plasma transfer set, 18" long with two spike couplers and a roller clamp	1TC*00505	100/case
Sampling site coupler with easy to pierce membrane for needle sampling	1TC*MP1	50/case



MEDISYSTEMS Apheresis Needles with MasterGuard® Anti-Stick Needle Protector

Product Description	Catalog No.	
16G x 1" needle, long bevel cannula with backeye, 12" tubing with clamp	8-P9-7006MGLB	250/case
17G x 1" needle, long bevel cannula with backeye, 12" tubing with clamp	8-P9-7007MGLB	250/case
17G x 1 1/4" needle, with backeye, 12" tubing without clamp	8-P9-7117MG	250/case
18G x 1 1/4" needle, with backeye, 12" tubing without clamp	8-P9-7118MG	250/case
18G x 1" needle, with backeye, 12" tubing with clamp	8-P9-7008MG	250/case



Blood Collection Accessories (Latex Free)

Product Description	Catalog No.			
VENOJECT® Multi-sample luer adapter, use with Standard Tube Holder and SampLok ® Tube Holder. 30 boxes containing 100 pieces each.	1XX*MN2000T	3000/case		
VENOJECT Standard Tube Holder, use with multi-sample luer adapter. For use with Diversion Blood Sampling Arm® and Blood Sampling Arm blood bags to collect samples. 20 bags containing 50 pieces each.	3P-1316R	1000/case		
DonorCare ® Needle Guard for use with blood collection sets. Safety device to reduce risk of needlestick injury.	8-A100100	1000/case		
SampLok Tube Holder for use with multi-sample luer adapter. For use with blood bags with Diversion Blood Sampling Arm and Blood Sampling Arm to collect samples. DonorCare can be locked in holder after phlebotomy.	8-A100400	1000/case		
SampLok Kit- Tube Holder with Luer Adapter. For use with Diversion Blood Sampling Arm and Blood Sampling Arm blood bags to collect samples. DonorCare can be locked in holder after phlebotomy. Sterile, individually packaged.	8-A100500	350/case		
Phlebotomy Needle Cap Holder. Holds needle cap while replacing on needle.	8P-7624A	50/case		
Multi-function Hand Stripper designed to strip tubing, crimp aluminum sealing clips and cut tubing.	Call for info	1/case		
Aluminum Sealing Clips to be used with the Multi-function Hand Stripper to permanently seal blood bag tubing.	1ME*ACS121	1000/case		
Plasma Separation Stand.	1ME*ACS201	1/case		
Green Tubing Clips. Plastic open end clip, slips on blood bag tubing to clamp close.	1BB*CLIPS	500/case		
Blue Slide Clamp. Plastic looped end clip, slips on blood bag tubing to clamp close.	1TF*03X001	500/case		
Soft, green hand squeezer.	8BB*SQUEEZER	25/case		
Small, multi-pocket organizer to hold blood collection accessories. Washable.	8BB*ORGANIZER	1/case		