



www.GenesisBPS.com
1-866-712-5663

Blood Collection & Processing Accessories



LABath Clean™

CAT# 402-05

- Water Bath Clarifier
- 240 mL concentrate



Plasma Guard Thawing Bags™

CAT# 401-03

- Protects ports during thawing process
- Unique double seal prevents/eliminates leaks



Sure Seal Thawing Bags™

CAT# 401-021

- Protects ports during thawing process
- Strong seal prevents/eliminates leaks



4" Blood Shoe™ Blood Bag Holders

CAT# 400-354-A

- Positions blood bags upright
- Labels are legible
- Clear, durable and easily cleaned



6" Blood Shoe™ Blood Bag Holders

CAT# 400-356-A

- Positions blood bags upright
- Labels are legible
- Clear, durable and easily cleaned



GENESIS Tube Stripper

CAT# 4A-HS-002

- High quality aluminum alloy
- Lightweight and sturdy
- Ergonomically designed

IV ADMIN SETS

CAT# 60150-A (Vented)

CAT# 60150C-A (Universal)

CAT# 60151C-A (Pediatric)

CAT# 60153C-A (Universal with 15 micron filter)

CAT# 60160-A (Non-Vented)

Blood Bags & Phlebotomy Supplies

GenesisBPS carries a full line of equipment and supplies to meet all your collection needs, including:

- Blood collection bags/transfer bags
- Sealing clips/grommets
- Hand/phlebotomy squeezers

GENESIS BPS™

GenesisBPS meets client needs worldwide for the collection, processing, storage and preparation of blood and other biological fluids. We do so by combining advanced technology, quality standards and operational efficiencies, working together to provide products and services of outstanding precision and reliability.

www.GenesisBPS.com
1-866-712-5663



GENESIS™

BLOOD PRODUCT SPECIALTIES

A Unique Line of Supplies for
Blood Bank & Cellular Therapies



GENESIS BPS™

Where the global bloodstream beginsSM



**Blood Bank Specialty
Equipment and Supplies**

GENESIS RAPIDWELD™ II STW sterile tube welder



CAT# 480-RW-2607-02

A technologically advanced, compact and waferless sterile tube welder

- Efficient & eco-friendly
- Waferless sterile welding
- Advanced blade module
- Touch screen interface
- Infrared heating element
- Patented tube holder clamps

Pending 510(k), not available for sale within the United States.



GENESIS™ DCM 3000™ BLOOD COLLECTION MIXER

CAT#425-CM760

Setting the standard for data collection

The DCM 3000 is an accurate and reliable mixer designed for fixed site and mobile collections.

- Blood Collection Data Management System
- A/C or rechargeable battery
- User-friendly interface
- ISBT 128 compatible barcode reader



GENESIS™ BLOOD COLLECTION MIXER

CAT# 425-CM735A

An accurate, durable and reliable automated mixer

- Weight accuracy +/- 1 gram
- Automatic taring
- Rechargeable battery



TCD® STERILE TUBE WELDER CAT# 480-B40-00-1245

- Small footprint
- Single wafer loading—no jamming
- No validation test kits required

Blood Bag Temperature Monitors



BLOOD TEMP IO

CAT# 3-BT10

- 10°C monitoring
- Monitor transportation



BLOOD TEMP 6

CAT# 3-BT06

- 6°C monitoring
- Monitor storage

Therapeutic Phlebotomy Products

MACOPHARMA

Therapeutic Phlebotomy Sets

CAT# 4M-VSL 7001PD CAT# 4M-VSL 7000YD

- Improve standard of care
- Graduated therapeutic phlebotomy bags



GENESIS™ CM350

Blood Collection Monitor

CAT# 425-CM350

- Electronically monitors and measures blood flow and volume
- Automatic weight calibration
- Automatically clamps with pre-set volume



Terumo Dry Collection Bag

CAT# 4T-1BB*D606A

- 600mL capacity
- 16 gauge ultra-thinwall donor needle



DonorCare Needle Guard

CAT# 4T-8-A100100

- For use with blood collection sets
- Reduces risk of needle stick injury



STEMIX MIXING AND COOLING

CAT# 43-CCS2017

- Temperature controlled mixing of cryoprotectant
- Compact design with minimal footprint
- Easy to operate, clean and maintain



GENESIS™ TUBE SEALERS

Innovative multi-head designs customized to fit your needs.

The SE530 & SE540 sealers can be used stand alone or linked as a multi-head sealing unit.

Multi-Head Sealer



RAPID SEAL II™ SE530

CAT# 428-SE530

- Space saving benchtop sealer
- Fixed sealing head



RAPID SEAL II™ SE540

CAT# 428-SE540

- Space saving benchtop sealer
- Fixed and hand-held sealing unit



RAPID SEAL II™ SE330

CAT# 428-SE330

- Standard benchtop sealer
- Fixed sealing head



RAPID SEAL II™ SE340

CAT# 428-SE340

- Multiple station benchtop sealer
- Fixed and hand-held sealing unit



RAPID SEAL II™ SE640

CAT# 428-SE640

- Battery powered portable sealing unit
- Hand-held sealing unit

ZipThaw® DRY THAWING SYSTEM



CAT# 47-ZT2020

- Consistent temperature control
- Small footprint, light-weight and portable
- Rapidly and completely thaws specimens
- Reduces the risk of any cross contamination
- 3-year full replacement warranty

Cord Blood & Cellular Laboratory Products



ELECTRONIC OPTICAL PLASMA EXPRESSOR

CAT# 422-ES315C

- Audible alarm
- Optical sensor detects red cells
- Automatically clamps, separating plasma from red cells



CORD BLOOD SEALER

CAT# 428-SE470RS

- Specially designed sealing head
- Designed to seal the Pall™ cord blood freezing-bag chamber



SEALING JIG

CAT# 428-UCBJIG

- Designed to work with the Pall™ cord blood freezing-bag
- Designed to accurately guide tube sealing head



GENESIS™ CB220

Centrifuge Balance

CAT# 419-CB220

- Electronic differential digital scale
- Displays weight differential of 2 objects

To order, call Genesis BPS at

1-866-712-5663

The Genesis RapidWeld™ STW and TCD® are protected by multiple patents.