

GENESIS™

BLOOD COLLECTION MONITOR

MODEL# CM350

CAT# 425-CM350



Keep a close watch on your donation – automatically.

- Electronically monitors and measures blood flow and volume
- Automatically clamps with pre-set volume
- Automatic weight calibration
- Rechargeable battery
- Dual function keypad
- Audible/visual alarms
- Mix alarm reminder
- Programmable default settings



GENESIS BPS™

Technology is in our Blood.™

GENESIS™

BLOOD COLLECTION MONITOR

MODEL# CM350
CAT# 425-CM350



SPECIFICATIONS

Dimensions	6.4"L x 4.8"H x 2.4"W
.....	162L x 60W x 120H mm
Collection volume	0 - 1199mL range
Weight display	0 - 1199gm range
Target volume range	100 - 600mL range
Elapsed time display	0:00 - 19:59 (min:sec)

ACCURACY

Weight tolerance	0 to +3 grams
Weight accuracy	+/- 1 gram
Time display accuracy	1 second

SPECIAL FEATURES

- Enhanced digital display
 - Pre-set collection volume by operator
 - Collected volume in mL
 - Total weight in grams
 - Elapsed time of donation
- Audible/visual alarms
 - Mix alarm reminder
 - End of collection
 - Slow bleed/low flow
 - Low battery
 - Wrong key
- Three indication LED's
- Electronically monitors and measures blood flow & volume
- Automatically clamps with pre-set volume
- Unique design for selecting specific programs suitable for collection needs
- Performs automatic weight calibration
- Dual function keypad
- Ability to mute audible alarms
- Automatic shut off after 30 minutes of idle time conserves battery power
- Rechargeable NiMH cell 7.2V/1.8ah battery

BENEFITS

- Reduces over or under bleeds (QNS units)
- Automatically tares collection bag
- Small, compact, designed for mobile use
- Mounts on any 1/2" pole
- Optional accessories available
- One year limited warranty

POWER REQUIREMENTS

Rechargeable NiMH cell, 7.2V/1.8ah
A/C current

WEIGHT

Collection monitor.....	1.50 lb
Collection monitor with pole mount and A/C adaptor....	2.25 lb
Shipping weight	4.40 lb
Package	One unit per carton

ACCESSORIES INCLUDED

Pole mount
A/C adaptor
500 gram standard weight

OPTIONAL ACCESSORIES

500 gram standard weight
Stand (pole and pedestal)
Aluminum carrying case

To order, call Genesis BPS at

1-866-712-5663



GENESIS BPS™

Technology is in our Blood.™